

bs-19425R

[Primary Antibody]

PRRC2C Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) Reactivity: (predicted: Human, Mouse, Rabbit, Pig, Sheep, Cow, Zebrafish, Dog, Horse, Rhesus monkey, Orangutan, Elephant) Predicted MW.: 317 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 23215	SWISS: Q9Y520	
Target: PRRC2C		
Immunogen: KLH conjugated synthetic peptide derived from human PRRC2C: 1-100/2896.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		