

**bs-0310Rs**

**[ Secondary Antibodies ]**

## **Rabbit Anti-Chicken IgG H&L whole serum**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### **DATASHEET**

**Host:** Rabbit

**Isotype:** IgG

**Reactivity:** Chicken

**Clonality:** Polyclonal

**Target:** Rabbit Anti-Chicken IgG H&L whole serum

**Purification:** Unpurified

**Storage:** Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.