

bs-0310Rs

[Secondary Antibodies]

Rabbit Anti-Chicken IgG H&L whole serum



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Rabbit Anti-Chicken IgG H&L whole serum

Purification: Unpurified

Storage: Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.

Applications: **WB** (1:10000-50000)

IHC-P (1:100-500)

IHC-F (1:10000-50000)

IF

ELISA (1:10000-100000)

Reactivity: Chicken