

bs-0373R

[Secondary Antibodies]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Rabbit Anti-Dog IgM

— DATASHEET —

Host: Rabbit

Isotype: IgG

Applications: **IHC-P** (1:100-500)

IHC-F (1:100-500)

IF (1:50-200)

ELISA (1:500-1000)

Clonality: Polyclonal

Target: Rabbit Anti-Dog IgM

Immunogen: Native Dog IgM: full length.

Purification: affinity purified by Protein A

Storage: 0.01M PBS (pH7.4).

Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.

Reactivity: Dog