

## AMH Mouse mAb

|                      |                                                                                 |
|----------------------|---------------------------------------------------------------------------------|
| Cat. No:             | V4902                                                                           |
| Name:                | AMH Mouse mAb                                                                   |
| Storage:             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:          | Monoclonal                                                                      |
| Clone No:            | 1D5                                                                             |
| Host/Source:         | Mouse                                                                           |
| Subtype:             | IgG2a                                                                           |
| Reactivity:          | (predicted:&nbsp;Human)                                                         |
| Format:              | Liquid                                                                          |
| Buffer:              | 0.01M PBS (pH7.4).                                                              |
| Applications:        | ELISA, ICA                                                                      |
| Antigen:             | Recombinant Human AMH Protein ( <a href="#">bs-43520P</a> )                     |
| Pairing Information: |                                                                                 |

| ELISA   |       | Detection |       |       |       |       |       |       |
|---------|-------|-----------|-------|-------|-------|-------|-------|-------|
|         |       | V4901     | V4902 | V4903 | V4904 | V4905 | V4906 | V4907 |
| Capture | V4901 | -         | +     | +++   | -     | -     | -     | -     |
|         | V4902 | -         | -     | +++   | -     | -     | -     | -     |
|         | V4903 | -         | -     | -     | -     | -     | -     | -     |
|         | V4904 | -         | -     | -     | -     | +++   | -     | -     |
|         | V4905 | -         | -     | -     | -     | -     | -     | -     |
|         | V4906 | -         | -     | -     | -     | -     | -     | -     |
|         | V4907 | -         | -     | -     | -     | -     | -     | -     |
| ICA     |       | Detection |       |       |       |       |       |       |
|         |       | V4901     | V4902 | V4903 | V4904 | V4905 | V4906 | V4907 |
| Capture | V4901 | -         | -     | -     | -     | -     | -     | -     |
|         | V4902 | -         | -     | -     | -     | -     | -     | -     |
|         | V4903 | -         | -     | -     | -     | -     | -     | -     |
|         | V4904 | -         | -     | -     | -     | -     | -     | -     |
|         | V4905 | -         | -     | -     | -     | -     | -     | -     |
|         | V4906 | +++       | -     | -     | -     | -     | -     | -     |
|         | V4907 | +         | +++   | -     | -     | -     | -     | -     |