
Phosphoserine, BSA Conjugated

Catalog Number: bs-11993PB

Concentration: 1mg/ml

AA Seq: Coupling protein

Activity: No

Form: Lyophilized or Liquid

Storage: 0.01M TBS(pH7.4).

Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Protein phosphorylation provides a signalling system that can be thought of as a kind of protein on/off switch for many cellular signalling pathways. Phosphorylation is observed on serine, threonine, tyrosine and histidine residues. Cellular networks underlying phosphorylation can be very complex and often occurs on multiple distinct sites on a given protein. Phospho-specific antibodies are becoming critical reagents both for basic research and for clinical diagnosis.