
Ulex Europaeus Lectin 1 Rabbit pAb

Catalog Number: bs-10025R

Target Protein: Ulex Europaeus Lectin 1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),
ELISA (1:5000-10000)

Reactivity: (predicted:Ulex europaeus)

Predicted MW: 27 kDa

Source: KLH conjugated synthetic peptide derived from Ulex Europaeus Lectin 1: 101-200/243.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

PRODUCT SPECIFIC PUBLICATIONS

[IF=10.6] Liu Jianyu. et al. New insights into allergic rhinitis treatment: MSC nanovesicles targeting dendritic cells. J NANOBIOTECHNOL.
2024 Dec;22(1):1-33 WB ; Mouse . 39294599