

**bs-10025R**

**[ Primary Antibody ]**

## Ulex Europaeus Lectin 1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> <b>WB</b> (1:500-2000) <b>IHC-P</b> (1:100-500) <b>IHC-F</b> (1:100-500) <b>IF</b> (1:100-500) <b>ICC/IF</b> (1:100-500) <b>ELISA</b> (1:5000-10000)  <b>Reactivity:</b> (predicted: Ulex europaeus)  <b>Predicted MW.:</b> 27 kDa
<b>Clonality:</b> Polyclonal		
<b>Target:</b> Ulex Europaeus Lectin 1		
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from Ulex Europaeus Lectin 1: 101-200/243.		
<b>Purification:</b> affinity purified by Protein A		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

### — SELECTED CITATIONS —

- **[IF=10.6]** Liu Jianyu. et al. New insights into allergic rhinitis treatment: MSC nanovesicles targeting dendritic cells. J NANOBIOTECHNOL. 2024 Dec;22(1):1-33 WB ;Mouse. 39294599