
Recombinant human EGF protein, GST

Catalog Number: bs-2049P

AA Seq: 971-1023/1207

Predicted MW: 32.2

Detected MW: 32 kDa

Tags: GST

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Storage: Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Epidermal growth factor (EGF) is an acid- and heat-stable 53 amino acid protein originally found in rodents and humans. It has been shown to be a potent mitogen for a variety of cell types both in vivo and in vitro. EGF binds to the EGF receptor on the surface of cells and mediates intrinsic phosphorylation of the receptor on tyrosine residues. It has been detected in nearly all body fluids, such as urine (urogastrone), saliva, milk and platelet-rich plasma. EGF, TGF β and vaccinia virus growth factor exhibit 30-40% amino acid homology. Several additional members of the EGF/TGF family have been described; these include Cripto, Amphiregulin and the heparin-binding EGF-like growth factor. Amphiregulin and the heparin-binding EGF-like growth factor both bind to the EGF receptor.

PRODUCT SPECIFIC PUBLICATIONS

[IF=8] Zhang Linhao. et al. Disruption of cholangiocyte-B cell crosstalk by blocking the CXCL12-CXCR4 axis alleviates liver fibrosis. CELL MOL LIFE SCI. 2023 Dec;80(12):1-19 Other ; . 38010435