
Recombinant HHV-5 UL16 protein, His

Catalog Number: bs-49094P

Concentration: >0.5 mg/ml

AA Seq: 22-184/230

Predicted MW: 18.1

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0) with 150mM NaCl.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: No data available.