

bs-6836R**[Primary Antibody]****GAS2L1 Rabbit pAb**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: Rat (predicted: Human, Mouse, Pig, Sheep, Cow, Dog, Horse) Predicted MW.: 73 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 10634	SWISS: Q99501	
Target: GAS2L1		
Immunogen: KLH conjugated synthetic peptide derived from Human GAS2L1: 201-300/681.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: The protein encoded by this gene, a member of the GAS2 family, is similar in sequence to the mouse protein Gas2, an actin-associated protein expressed at high levels in growth-arrested cells. Expression of the mouse Gas2 gene is negatively regulated by serum and growth factors. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2008].		