

## Ovalbumin, Cy3 conjugated

Catalog Number: bs-0283P-Cy3

Concentration: 1mg/ml

Form: Liquid


Storage: 10 mM TBS (pH=7.4) with 0.02% Proclin300 and 50% glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Background:** Chicken Ovalbumin is the major protein in the "white" of the egg (and a favorite antigen in immunological research). Egg white contains a variety of proteins including ovalbumin, conalbumin, ovomucoid and lysozyme. It belongs to the serpin family and the Ov serpin subfamily. Ovalbumin can cause an allergic reaction in humans.

Ovalbumin has been implicated in the development of the egg shell. Immunohistochemistry revealed that ovalbumin is found only in the mammary bodies of decalcified shell, and is not distributed throughout the shell matrix. These results indicate that ovalbumin is present during the initial phase of shell formation and becomes incorporated into the protein matrix of the mammary bodies. However, it is not yet clear whether ovalbumin at this site plays a specific role in shell mineralisation.

### VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.