

---

## PAR-2 OP Peptide

Catalog Number: bs-0436P

AA Seq: FSLLRY-NH<sub>2</sub>

Activity: Yes

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Selective PAR2 peptide antagonist. Reverses taxol-induced mechanical allodynia, heat hyperalgesia and PKC activation in ICR mice. Blocks ERK activation and collagen production in isolated cardiac fibroblasts. Also reduces symptoms in a mouse model of dermatophyte-associated itch.