
Gentamicin Rabbit pAb

Catalog Number: bs-0972R

Target Protein: Gentamicin

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:Gentamicin)

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Gentamicin is an aminoglycoside antibiotic, used to treat many types of bacterial infections, particularly those caused by Gram-negative bacteria. However, gentamicin is not used for *Neisseria gonorrhoeae*, *Neisseria meningitidis* or *Legionella pneumophila* bacterial infections (because of the risk of the patient going into shock from lipid A endotoxin found in certain gram negative organisms).