

bs-16139R

[Primary Antibody]

FLJ40504 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 284085 Target: FLJ40504 Immunogen: KLH conjugated synthetic peptide derived from human FLJ40504: 101-200/259. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: Highly similar to Keratin, type I cytoskeletal 18.	Isotype: IgG SWISS: Q8N7Q0 Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human) Predicted MW.: 26 kDa
---	---