

AMH Mouse mAb

Cat. No: V4905
Name: AMH Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 3E11
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA
Antigen: Recombinant Human AMH Protein ([bs-43520P](#))
Pairing Information:

ELISA		Detection						
		V4901	V4902	V4903	V4904	V4905	V4906	V4907
Capture	V4901	-	+	+++	-	-	-	-
	V4902	-	-	+++	-	-	-	-
	V4903	-	-	-	-	-	-	-
	V4904	-	-	-	-	+++	-	-
	V4905	-	-	-	-	-	-	-
	V4906	-	-	-	-	-	-	-
	V4907	-	-	-	-	-	-	-
ICA		Detection						
		V4901	V4902	V4903	V4904	V4905	V4906	V4907
Capture	V4901	-	-	-	-	-	-	-
	V4902	-	-	-	-	-	-	-
	V4903	-	-	-	-	-	-	-
	V4904	-	-	-	-	-	-	-
	V4905	-	-	-	-	-	-	-
	V4906	+++	-	-	-	-	-	-
	V4907	+	+++	-	-	-	-	-