

AMH Mouse mAb

Cat. No:	V4907
Name:	AMH Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6A9
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ICA
Antigen:	KLH conjugated synthetic peptide derived from human AMH (bs-43520P)
Pairing Information:	

ELISA		Detection						
		V4901	V4902	V4903	V4904	V4905	V4906	V4907
Capture	V4901	-	+	+++	-	-	-	-
	V4902	-	-	+++	-	-	-	-
	V4903	-	-	-	-	-	-	-
	V4904	-	-	-	-	+++	-	-
	V4905	-	-	-	-	-	-	-
	V4906	-	-	-	-	-	-	-
	V4907	-	-	-	-	-	-	-
ICA		Detection						
		V4901	V4902	V4903	V4904	V4905	V4906	V4907
Capture	V4901	-	-	-	-	-	-	-
	V4902	-	-	-	-	-	-	-
	V4903	-	-	-	-	-	-	-
	V4904	-	-	-	-	-	-	-
	V4905	-	-	-	-	-	-	-
	V4906	+++	-	-	-	-	-	-
	V4907	+	+++	-	-	-	-	-