

Alpha 1 microglobulin Mouse mAb

Cat. No: V6101
Name: Alpha 1 microglobulin Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 4A9
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant Human Alpha 1 Microglobulin Protein ([bs-43007P](#))
Pairing Information:

ELISA		Detection				
		V6101	V6102	V6105	V6106	V6107
Capture	V6101	-	-	-	-	-
	V6102	+++	-	-	-	-
	V6105	-	-	-	-	-
	V6106	+++	-	-	-	-
	V6107	-	-	-	-	-
ICA		Detection				
		V6101	V6102	V6105	V6106	V6107
Capture	V6101	-	-	-	-	+++
	V6102	-	-	-	-	-
	V6105	+++	-	-	-	-
	V6106	+	-	-	-	-
	V6107	-	-	-	-	-