
Alpha 1 microglobulin Mouse mAb

Cat. No:	V6104
Name:	Alpha 1 microglobulin Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4H3
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human Alpha 1 Microglobulin Protein: (bs-43007P)