

---

## Recombinant dog gastrin-releasing peptide protein, C-hFc (HEK293)

Catalog Number: bs-43120P

Concentration: 1mg/ml

AA Seq: 24-147/147

Predicted MW: 40.4

Tags: C-hFc

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: PBS (pH=7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.