

beta 2 Microglobulin Mouse mAb

Cat. No:	V0706
Name:	beta 2 Microglobulin Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4C8
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant human Beta-2-microglobulin (bs-43094P)
Pairing Information:	

ELISA		Detection			
		V0703	V0705	V0706	V0707
Capture	V0703	-	-	+++	-
	V0705	+++	-	-	-
	V0706	-	-	-	-
	V0707	+	-	-	-
ICA		Detection			
		V0703	V0705	V0706	V0707
Capture	V0703	-	-	-	+++
	V0705	-	-	-	-
	V0706	-	-	-	-
	V0707	+++	-	-	-