

beta 2 Microglobulin Mouse mAb

Cat. No: V0707
Name: beta 2 Microglobulin Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 4C11
Host/Source: Mouse
Subtype: IgG2b
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant human Beta-2-microglobulin ([bs-43094P](#))
Pairing Information:

ELISA		Detection			
		V0703	V0705	V0706	V0707
Capture	V0703	-	-	+++	-
	V0705	+++	-	-	-
	V0706	-	-	-	-
	V0707	+	-	-	-
ICA		Detection			
		V0703	V0705	V0706	V0707
Capture	V0703	-	-	-	+++
	V0705	-	-	-	-
	V0706	-	-	-	-
	V0707	+++	-	-	-