
SGTA Rabbit pAb

Catalog Number: bs-7830R

Target Protein: SGTA

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat)

Predicted MW: 34 kDa

Subcellular Cytoplasm ,Nucleus

Locations:

Entrez Gene: 6449

Swiss Prot: O43765

Source: KLH conjugated synthetic peptide derived from human SGTA: 181-280/313.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Co-chaperone that binds directly to HSC70 and HSP70 and regulates their ATPase activity.