
Recombinant mStrawberry protein, His

Catalog Number: bs-33175P

Concentration: >1 mg/ml

AA Seq: 1-236/236

Predicted MW: 26

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH=8.0) with 150mM NaCl

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.



Recombinant mStrawberry protein