
Recombinant human HSA protein (Active, HEK293)

Catalog Number: bs-48030P

Concentration: >0.5 mg/ml

AA Seq: 1-609/609

Detected MW: 60-65 kDa

Tags: Tag free

Activity: Yes

Endotoxin: <0.5 EU/μg as determined by LAL.

Purity: >95% as determined by SDS-PAGE

Purification: AC

Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. In general: 12 months from date of receipt, -20 to -80°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 3 months, -20 to -80°C under sterile conditions after reconstitution.

Background: Recombinant Human Serum Albumin is produced in engineered human cells cultured in animal product and serum free chemically defined media. HSA has ultralow endotoxin and is ideal for use as a carrier protein in place of BSA.