
Recombinant human MYOD1 protein, C-His

Catalog Number: bs-42036P

Concentration: >0.5mg/ml

Species: Human

AA Seq: 1-320/320

Predicted MW: 38

Tags: C-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a nuclear protein that belongs to the basic helix-loop-helix family of transcription factors and the myogenic factors subfamily. It regulates muscle cell differentiation by inducing cell cycle arrest, a prerequisite for myogenic initiation. The protein is also involved in muscle regeneration. It activates its own transcription which may stabilize commitment to myogenesis. [provided by RefSeq, Jul 2008]

VALIDATION IMAGES



The purity of the protein is greater than 40% as determined by reducing SDS-PAGE.