
NRTN Antibody Blocking Peptide

Catalog Number: bs-0073P

AA Seq: human Neurturin (151-197aa)

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Human Neurturin (NTN) and Human GDNF comprise a family of TGF-beta related neurotrophic factor which have trophic influences on a variety of neuronal populations. Neurturin promotes the survival of certain sympathetic and sensory neurins through interaction with distinct set of GDNF-like receptors