
Camk1g Antibody Blocking Peptide

Catalog Number: bs-1809P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: CaMKI gamma is a Calcium/calmodulin-dependent protein kinase which belongs to a proposed calcium-triggered signaling cascade. In vitro it phosphorylates transcription factor CREB1. It has been proposed that CaMKI gamma is involved in dendritic outgrowth mediated by CaM-dependant protein kinase kinase.