

CRP Mouse mAb

Cat. No: V0101
Name: CRP Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 2D5
Host/Source: Mouse
Subtype: IgG2a
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Human CRP protein ([bs-41211P](#))
Pairing Information:

ELISA		Detection						
		V0101	V0102	V0103	V0106	V0107	V0109	V0110
Capture	V0101	-	+++	-	-	-	-	-
	V0102	-	-	-	-	-	-	-
	V0103	-	-	-	-	-	-	-
	V0106	+++	-	-	-	-	-	-
	V0107	-	-	-	-	-	-	-
	V0109	-	-	-	-	-	-	-
	V0110	-	-	-	-	-	-	-
ICA		Detection						
		V0101	V0102	V0103	V0106	V0107	V0109	V0110
Capture	V0101	-	-	-	-	-	-	-
	V0102	+	-	-	-	-	-	-
	V0103	-	-	-	+++	-	-	-
	V0106	-	-	-	-	-	-	-
	V0107	-	-	-	-	-	+++	-
	V0109	-	-	-	-	-	-	-
	V0110	-	-	-	-	+	-	-