

CRP Mouse mAb

Cat. No:	V0102
Name:	CRP Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5D12
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Human CRP protein (bs-41211P)
Pairing Information:	

ELISA		Detection						
		V0101	V0102	V0103	V0106	V0107	V0109	V0110
Capture	V0101	-	+++	-	-	-	-	-
	V0102	-	-	-	-	-	-	-
	V0103	-	-	-	-	-	-	-
	V0106	+++	-	-	-	-	-	-
	V0107	-	-	-	-	-	-	-
	V0109	-	-	-	-	-	-	-
	V0110	-	-	-	-	-	-	-
ICA		Detection						
		V0101	V0102	V0103	V0106	V0107	V0109	V0110
Capture	V0101	-	-	-	-	-	-	-
	V0102	+	-	-	-	-	-	-
	V0103	-	-	-	+++	-	-	-
	V0106	-	-	-	-	-	-	-
	V0107	-	-	-	-	-	+++	-
	V0109	-	-	-	-	-	-	-
	V0110	-	-	-	-	+	-	-