

bs-6513R

[Primary Antibody]

A2BP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 54715</p> <p>Target: A2BP1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human A2BP1/Fox1: 151-250/397.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: RNA-binding protein that regulates alternative splicing events by binding to 5'-UGCAUGU-3' elements. Regulates alternative splicing of tissue-specific exons and of differentially spliced exons during erythropoiesis. Predominantly expressed in muscle and brain.</p>	<p>Isotype: IgG</p> <p>SWISS: Q9NWB1</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Cow, Dog, Horse)</p> <p>Predicted MW.: 44 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p>
---	--	---