
AD7c-NTP Mouse mAb

Cat. No:	V3006
Name:	AD7c-NTP Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2B4
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human AD7c-NTP Protein (bs-41194P)