

---

## Recombinant HHV-8 K3 protein

Catalog Number: bs-0864P

Concentration: >0.5 mg/ml

AA Seq: Unknown

Predicted MW: 37

Detected MW: 37 kDa

Tags: Tag free

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized

Storage: Lyophilized from 0.22um filtered solution in PBS (pH7.4) with 5mM DTT. Normally 5% trehalose is added as protectant before Lyophilization.

Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.