
Recombinant CMV L-CM protein

Catalog Number: bs-49064P

Concentration: >0.5 mg/ml

AA Seq: 1-110/110

Predicted MW: 25

Detected MW: 25 kDa

Tags: Tag free

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.