
Recombinant human PCSK9 protein, C-His (HEK293)

Catalog Number: bs-43105P

Concentration: >0.5 mg/ml

AA Seq: 31-692/692

Predicted MW: 73

Tags: C-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: May be implicated in the differentiation of cortical neurons and may play a role in cholesterol homeostasis.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.