

Recombinant mCherry protein, His

Catalog Number: bs-33131P

Concentration: >0.5 mg/ml

AA Seq: 1-236/236

Predicted MW: 25

Detected MW: 33 kDa

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: mCherry is a fluorophore (a fluorescent molecule) used in biotechnology as a tracer to follow the flow of fluids, as a marker when tagged to molecules and cells components. mCherry is a monomeric fluorescent construct with peak absorption/emission at 587 nm and 610 nm, respectively. It is resistant to photobleaching and is stable. mCherry is sometimes preferred to other fluorophores due to its colour, as well as its photostability compared to other monomeric fluorophores.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.