

---

## Recombinant mOrange protein, His

Catalog Number: bs-33183P

Concentration: >1 mg/ml

AA Seq: 1-236/236

Predicted MW: 26

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-Hcl (pH=8.0) with 150 NaCl

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

### VALIDATION IMAGES

---



The purity of the protein is greater than 85% as determined by reducing SDS-PAGE.



Recombinant mOrange protein