
ANGPTL2 Antibody Blocking Peptide

Catalog Number: bs-0537P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ANGPTL2/ANG-2 (Angiopoietin-like Protein 2) Induces sprouting in endothelial cells through an autocrine and paracrine action. [Subcellular location] Secreted. [Tissue specificity] Widely expressed in heart, small intestine, spleen and stomach. Also found in lower levels in colon, ovary, adrenal gland, skeletal muscle and in prostate. [PTM] N-glycosylated. [Similarity] Contains 1 fibrinogen C-terminal domain. ANGPTL2/ANG-2 (Angiopoietin-like Protein 2).