

---

## CCL22 Antibody Blocking Peptide

Catalog Number: bs-1761P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: CCL22; DC/B-CK; A-152E5.1; ABCD 1; ABCD-1; ABCD1; CC chemokine STCP-1; ccl 22; MGC34554; SCYA22; small inducible cytokine subfamily A (Cys-Cys) member 22; STCP 1; STCP-1; MDC; stimulated T cell chemotactic protein 1; CCL22\_Mouse; C-C motif chemokine 22; MDC(1-69); Macrophage-derived chemokine; Small-inducible cytokine A22; Stimulated T-cell chemotactic protein 1; MDC(3-69); MDC(5-69); MDC(7-69); CCL22\_HUMAN.