
Recombinant human IgG2 Fc protein,(HEK293)

Catalog Number: bs-43576P

Concentration: >1mg/ml

Species: Human

AA Seq: 99-326/326

Predicted MW: 25.7

Tags: Tag Free

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH=7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.