

---

## Recombinant human Ferritin Light Chain, His

Catalog Number: bs-41189P

Concentration: >0.5 mg/ml

AA Seq: 1-175/175

Predicted MW: 19

Detected MW: 27 kDa

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid


Storage: 20mM Tris-HCl (pH8.0) with 150mM NaCl.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** This gene encodes the light subunit of the ferritin protein. Ferritin is the major intracellular iron storage protein in prokaryotes and eukaryotes. It is composed of 24 subunits of the heavy and light ferritin chains. Variation in ferritin subunit composition may affect the rates of iron uptake and release in different tissues. A major function of ferritin is the storage of iron in a soluble and nontoxic state. Defects in this light chain ferritin gene are associated with several neurodegenerative diseases and hyperferritinemia-cataract syndrome. This gene has multiple pseudogenes. [provided by RefSeq, Jul 2008].

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.