
Recombinant human RCBTB2 protein, His

Catalog Number: bs-49038P

Concentration: >0.5 mg/ml

AA Seq: 1-551/551

Predicted MW: 60

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0) with 150mM NaCl.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a protein containing two C-terminal BTB/POZ domains that is related to regulator of chromosome condensation (RCC). The encoded protein may act as a guanine nucleotide exchange factor. This gene is observed to be lost or underexpressed in prostate cancers. There is a pseudogene of this gene on chromosome 10. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2013]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.