
penicillin/BSA

Catalog Number: bs-10687PB

Concentration: 1mg/ml

AA Seq: Coupling protein

Tags: N/A

Activity: No

Endotoxin: Not analyzed

Form: Liquid

Storage: 0.01M TBS(pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Penicillin is a group of Beta-lactam antibiotics used in the treatment of bacterial infections caused by susceptible, usually Gram-positive, organisms. β -lactam antibiotics work by inhibiting the formation of peptidoglycan cross-links in the bacterial cell wall, which results in cytolysis.