

bs-0245R

[Primary Antibody]

ELP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: ELP1

Immunogen: KLH conjugated synthetic peptide derived from Arabidopsis thaliana IKI3 family protein: 1256-1319/1319.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

Reactivity: (predicted: Arabidopsis
Thaliana)

**Predicted
MW.:** 147 kDa