
Recombinant human AD7c-NTP protein, GST

Catalog Number: bs-41187P

Concentration: >0.5 mg/ml

AA Seq: 1-200/375

Predicted MW: 22

Detected MW: 45 kDa

Tags: GST

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC


Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0) with 150mM NaCl.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: AD7c-NTP is detected in increased concentration in the cortical neurons, brain-tissue extracts, cerebrospinal fluid, and urine early in the course of AD neurodegeneration, and its level is positively correlated with the severity of dementia. All these characteristics make it a possible biomarker for AD.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.