

Recombinant chicken MIER3L protein, N-His

Catalog Number: bs-49157P

Concentration: >0.5 mg/ml

AA Seq: 264-419/640

Predicted MW: 18.3

Detected MW: 25 kDa

Tags: N-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE


Purification: AC

Form: Liquid

Storage: 2M Urea with 20mM Tris-HCl and 15mM L-Arginine.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.