
Recombinant human Mesothelin, N-His

Catalog Number: bs-42016P

Concentration: >0.5 mg/ml

AA Seq: 296-580/630

Predicted MW: 32.4

Tags: N-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 4M Urea (pH9.0) with 20mM Tris-HCl.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Mesothelin(MSLN) may play a role in cell adhesion and cell shape dynamics; human homolog highly expressed in renal cell carcinoma. malignant tumors, such as a skin lump, ovary cancer and the pancreas cancer...etc..Mainly used for the research of the malignant tumor.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.