

---

## Recombinant human EGFR Fragment protein, C-His

Catalog Number: bs-10695P

Concentration: >0.5 mg/ml

AA Seq: 178-304/1210

Predicted MW: 13

Detected MW: 15-17 kDa

Tags: C-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0) with 8M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor. Binding of the protein to a ligand induces receptor dimerization and tyrosine autophosphorylation and leads to cell proliferation. Mutations in this gene are associated with lung cancer. Multiple alternatively spliced transcript variants that encode different protein isoforms have been found for this gene. [provided by RefSeq, Jul 2010]