

---

## Recombinant human ICAM1 / CD54 protein, C-His

Catalog Number: bs-10697P

Concentration: >0.5 mg/ml

AA Seq: 321-420/532

Predicted MW: 13.6

Detected MW: 14 kDa

Tags: C-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 10mM TBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor. [provided by RefSeq, Jul 2008].

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.