
prox1 Antibody Blocking Peptide

Catalog Number: bs-2774P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Prox1 is a prospero-related homeobox gene. Its expression pattern suggests it has a role in a variety of embryonic tissues and may play a fundamental role in early development of the CNS. It may regulate gene expression and the development of postmitotic undifferentiated young neurons. Subcellular Location : Nucleus.